

Technical Data Sheet

Article-No.: **09283**

Description:

urine bottle
for men



DIMENSIONS

length	width	Ø cap
310 mm	85 mm	58 mm

PRODUCT DESCRIPTION

material bottle	☐ Polypropylene (PP)	☐ Polyethylene (PE)	☒ Polycarbonate (PC)	☐ Polyvinylchloride (PVC)			
material lid	☒ Polypropylene (PP)	☐ Polyethylene (PE)	☐ Polycarbonate (PC)				
colour bottle	☐ milky	☒ transparent					
colour lid	☐ blue	☒ white					
characteristics	☒ disinfectable	☒ autoclavable up to 130°C	☐ autoclavable up to 90°C	☒ angular	☐ autoclavable up to 60°C		

REGULATORY AFFAIRS

PPE-Regulation (EU) 2016/425	☐ Category I	☐ Category II	☐ Category III	☒ no PPE-article			
MD-Regulation (EU) 2017/745	☐ Class I	☐ Class II	☐ Class III	☐ sterile	☐ measuring function	☒ no medical device	

LOGISTIC DATA SUBPACKING

general information		
material	bag	
pieces per subpacking	1	
EAN subpacking	-	
measures & size		
length	450 mm	
width	180 mm	
height	-	
weights		
size	net weight	brut weight
uni	140 g	145 g

LOGISTIC DATA PALETTE

general information		
kind of palette		euro-palette
measures & size		
cartons per layer		4
layer per palette		5
height of the palette		225 cm
weights		
size	net weight	brut weight
uni	83 kg	108 kg

LOGISTIC DATA OUTER PACKING

general information	
material	carton
subpackings per outer packing	25
EAN outer packing	4044941042206
measures & size	
length	500 mm
width	320 mm
heigh	420 mm
weights	
size	net weight
uni	3.625 g
	brut weight
	4.125 g

Technical Data Sheet

Article-No.: **09283**

Description:
**urine bottle
for men**



WARNINGS AND SAFETY INFORMATION

storage	Store items in original packaging in a cool and dry place without additional weight, protect from direct sunlight.
use and control	Check the product before first use and regularly for cracks or other damage and replace it if necessary. Close the lid properly after use to minimize odours. Do not dispose of the contents of the bottle in the environment. If necessary, transport the bottle in a sealed container. Wear gloves when handling the bottle and its contents.
disposal	Used bottles should be disposed of in accordance with the disposal regulations of the local waste disposal company. Unused products can be disposed of with household waste.
warnings / safety information	The urine bottle is intended exclusively for holding urine. Any other use can lead to hygienic or mechanical problems. Urine bottles are not designed to be used with other liquids or chemicals.

rev-no. 1
date 01.04.2025

changes and errors excepted