

Technical Data Sheet

Article-No.: L-09152

Description: Dispenser for 90 medicine cups



SIZES		
length	width	height
395 mm	50 mm	-

PRODUCT DESCRIPTION							
material	<input type="checkbox"/> stainless steel	<input type="checkbox"/> acrylic glass	<input type="checkbox"/> PVC coated wire	<input type="checkbox"/> ABS plastic	<input type="checkbox"/> casing: aluminium	<input type="checkbox"/> pump: plastic	<input checked="" type="checkbox"/> polymethylmethacrylat (PMMA)
colour	<input type="checkbox"/> silver	<input checked="" type="checkbox"/> transparent	<input type="checkbox"/> white	<input type="checkbox"/> blue			
characteristics	<input type="checkbox"/> with movable hooks	<input checked="" type="checkbox"/> with fixing material	<input type="checkbox"/> adjustable stop for depth adjustment		<input checked="" type="checkbox"/> for wall mounting		
desription							
suitable for	09174, 09176						

LOGISTIC DATA SUBPACKING		
general information		
material	#NV	
pieces per subpacking	0	
GTIN subpacking	-	
measures & sizes		
length	0 mm	
width	0 mm	
height	0 mm	
weights		
size	net weight	brut weight
uni	400 g	435 g

LOGISTIC DATA OUTER PACKING		
general information		
material	carton	
subpackings per outer packing	1	
GTIN outer packing	-	
measures & sizes		
length	420 mm	
width	90 mm	
height	90 mm	
weights		
size	net weight	brut weight
uni	435 g	935 g

LOGISTIC DATA PALLETTE		
general information		
kind of palette	euro-palette	
measures & sizes		
cartons per layer	-	
cartons per palette	300	
palette height	200 cm	
weights		
size	net weight	brut weight
uni	281 kg	306 kg

Technical Data Sheet

Article-No.: **L-09152**

Description: **Dispenser for 90 medicine cups**



WARNINGS AND SAFETY INFORMATION

storage	Store items in original packaging in a cool and dry place without additional weight, protect from direct sunlight.
warnings / safety information	<ul style="list-style-type: none">- Mount the product firmly and securely on the wall. Check regularly that the bracket is stable.- When installing the dispenser, ensure that the walls provide the necessary stability- Before use, ensure that all parts are securely mounted and intact.
cleaning instructions	Can be cleaned easily. Use suitable cleaning agents and soft cloths to avoid scratches.

rev-no.: 1
date: 24.03.2025

changes and errors excepted